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MPV TRUMA

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MicroDrop® Pro

Kompressor / Compressor

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
.....

Warning:

Remove the compressor from carrying case before use!

Components and Accessories - (Pic. A)

1. Main unit
2. ON / OFF switch
3. Power cord
4. Air outlet
5. Filter holder
6. Nebulizer holder
7. Air vent openings
8. Handle

Standard Accessories

9. *Pro* nebulizer
10. Mouth piece
11. Nose piece
12. Air tube
13. Adult mask
14. Child mask
15. Carrying case

Important:

**Read carefully instructions manual
before using this device.**

General Information

MicroDrop® *Pro* is a compressor nebulizer system for aerosol therapy for intensive use. MicroDrop® compressor is supplied with the *Pro* nebulizer, a new and patented exhalation/inhalation valves system, which increases the amount of inhaled medication, minimizing the medication waste.

Innovative Double-Valve System Operating Detail

Illustration B

During inhalation phase, the inhalation valve (pic. B / 16) opens, automatically-operated by patient's breath, thereby allowing an additional air flow (pic. B / 17) to the compressor-generated one (pic. B / 18), therefore increasing the nebulization of medication and the medication quantity effectively inhaled.

Illustration C

During exhalation phase, the inhalation valve (pic. B / 16) closes avoiding medication to be dispersed in the environment while the exhalation valve (pic. C / 19) opens thereby allowing the exhalation.

Precautions

1. **Only use the device for its intended purpose as described in this instruction manual, therefore as a compressor nebulizer system for aerosoltherapy, following the indications of your doctor.**

Any other form of use constitutes an improper use and is therefore dangerous. The manufacturer cannot be held liable for any damage caused by improper, incorrect or unreasonable use, or if the unit is connected to electrical installations which do not conform with safety regulations.

2. Keep this manual for future reference.
3. Do not operate the unit in presence of any anesthetic, inflammable mixtures or oxygen.
4. The correct device functioning can be affected by electromagnetic interference which exceeds the limits indicated by the European standards in force.
5. In case of failure and / or malfunction, read the "POSSIBLE PROBLEMS AND HOW TO SOLVE THEM" section. Do not handle nor open the compressor housing.
6. Any repairs must be carried out by an authorized dealer using only original spare parts. Safety of the device can be compromised in case the above mentioned indications are not met.
7. When using any electrical appliance certain important safety measures must always be observed, including the following:

- use only original accessories and components;
 - never submerge the unit in water;
 - the device is not protected against water penetration;
 - never touch the unit when your hands are wet or moist;
 - do not leave the unit outdoor;
 - the unit, when operating, must be placed on a stable and horizontal surface;
 - make sure the air vent openings are not obstructed;
 - do not allow unsupervised children or infirm persons to use the unit;
 - do not disconnect the unit by simply pulling out the plug from the wall socket.
8. Make sure that the electrical rating shown on the rating plate on the bottom of the unit corresponds to your main voltage and frequency before plugging in the device.
 9. If the power plug provided with the device does not fit your wall socket, have the plug changed by an electrical appliance dealer. Do not use any adapter, simple or multiple, and / or extension cable. In case their use is necessary, make sure they are in compliance with safety standards, paying attention they do not exceed the maximum load.
 10. Do not leave the unit plugged when not in use: unplug the device from the wall socket when not operating.
 11. Follow the manufacturer's instructions for installing the device. The manufacturer cannot be held responsible for any damage caused to persons, animals or things by incorrect installation.
 12. The power supply cord cannot be replaced by the user. In case the power supply cord becomes damaged, address to the technical servicing authorized by the manufacturer for replacement.
 13. The power supply cord should always be fully unwound in order to prevent dangerous overheating.
 14. Before performing any maintenance or cleaning operation, disconnect the device from the main supply, by unplugging it or by turning the main supply switch off.
 15. Some parts of the unit are so small that they could be swallowed by children; do not let unattended children have access to the unit.
 16. If you decide not to use the device any longer, it is recommended to dispose of it according to the regulations in force.
 17. Remember to:
 - use this device only with medications prescribed by your doctor;
 - make the treatment only using the accessory recommended by your doctor for the therapy.
 - use the “nose piece” accessory only if expressly indicated by your doctor and paying attention **NEVER** to introduce inside the nose the nasal bifurcation, but only bring it as close as possible.

Device Safety Systems

Thermal Protector

The inhalation device is provided with a temperature protection device which is released when the device is overheating switching the unit off.

Operating time:

90 min. ON, 15 min. OFF

Instruction for Use

Remove the compressor from carrying case before use.

Check the unit before each use, in order to check for possible function anomalies and / or damages due to transport and / or storage.

When inhaling, sit upright and relaxed at a table and not in armchair as otherwise you will constrict your respiratory passages and impair the effect of inhalation.

1. As soon as you have unpacked the device, check it for visible damages or defects. In particular, check for cracks in the plastic housing, which could expose electrical components.
2. Check for accessories integrity.
3. **Before using the device**, proceed with the cleaning operations as described in the "CLEANING AND DISINFECTION" section.
4. Open the nebulizer (pic. A / 9) by turning the top in an anticlockwise direction (pic. D). On the top an arrow, together with "CLOSE" sign, suggest how to close the nebulizer.
5. Make sure the water conduction cone is properly fitted onto the air conduction cone inside the nebulizer (pic. E).
6. Put the prescribed medication dose into the nebulizer (pic. F). The graduated scale can be useful to control the correct medication quantity.
7. Close the nebulizer by turning the two parts in a clockwise direction, paying attention that they are thoroughly sealed.
8. Connect one end of the air tube to the nebulizer (pic. G) and the other end to the air outlet on the device (pic. A / 4).
9. Use the *Pro* nebulizer with inhalation valve only with mouthpiece having exhalation valve. When using the masks or the nosepiece, it is not necessary to keep inhalation valve: remove it (pic. B / 16) and connect the nosepiece or the masks directly onto the nebulizer (pic. I).
10. Plug the device to the wall socket, making sure that the electrical rating shown on the rating plate on the bottom of the unit corresponds to your main voltage.
11. To start the treatment, switch the device on by pressing the switch (pic. A / 2) in the front of the device onto the "I" position.

12. Inhale the aerosol solution using the prescribed accessory.
13. When the treatment has been completed, switch the unit off by setting the switch position "O" and disconnect the plug from the wall socket.
14. Wash the nebulizer and accessories as described in the "CLEANING AND DISINFECTION OF THE ACCESSORIES" section.
15. Put the nebulizer, accessories and power cord inside the carrying case (pic. A / 15).

Cleaning and Maintenance

Cleaning of the main unit

At home:

Use a clean dry cloth and mild detergents for cleaning the unit.

Dealer:

After cleaning the device can be complemented for use with a new patient.

To do so, the following steps must be performed:

- wipe disinfection (with common disinfectants)
- add a new nebulizer set (e.g. the practical year set)
- pressure control of the compressor

The device is now ready for use by a new patient.

During cleaning operations, make sure the internal parts of the unit are not in contact with liquids and the power cord is not connected.

Cleaning and disinfection of the accessories

See nebulizer manual for detailed information.

After each treatment:

The nebulizer, the mouth piece and the nose piece must be sterilized by boiling in water (for max. 10 minutes) or in autoclave (max. 121°C). We suggest you disinfect the mask by using cold disinfecting liquids (solutions with hypochlorite). Dry the nebulizer using a clean cloth; after connect it to the air outlet (pic. A / 4), switch the device on and let it work freely for 5-10 minutes.

Do not boil nor autoclave the air tube and the masks.

Following carefully the cleaning and disinfecting instructions is very important for the product performance and for a successful therapy.

Replacement of the nebulizer

The nebulizer must be replaced after a long period of inactivity, or when it is obstructed by dry medication, dust, etc. Your nebulizer should last from 6 months to 1 year with general use. Use original nebulizers only. (Year Set *Pro* nebulizer M 53710-04)

Replacement of the filter

The MicroDrop® *Pro* Porex filter must be replaced approximately after 500 hours of use or each year; control the filter periodically (10 - 12 aerosol treatments); when it turns grey or brown or wet, replace it; open the filter holder (pic. A / 5) and replace the filter with a new one. Don't try to clean the filter and to reuse. Use original filters only.

Do not use the device without filter.

Possible Problems and how to solve them

If the device does not switch on:

1. make sure the plug is firmly fitted to the wall socket;
2. the thermal protector might have been released. The thermal protector may be released for the following reasons:
 - the device has been in operation for a too long time;
 - the device was used near heaters or in hot environments;

- the air vent openings (pic. A / 7) are obstructed and do not allow for the motor ventilation.

In this case, set the ON/OFF switch to "O" and let it rest for at least 15 minutes.

If the device does not nebulize, or only very weakly, check the following points:

1. make sure that the ends of the air tube are fitted tightly onto the main unit and onto the nebulizer;
2. check whether the nebulizer is empty from medication or full of an excessive medication quantity (max. 15 ml);
3. check whether the nebulizer nozzle is obstructed;
4. check whether the air conduction cone is inserted (pic. E).

Should the device not start working properly again, address to an authorized technical service centre.

Maintenance and Repairs

In case of failure, address to qualified personnel authorized by MPV TRUMA.

Never open the device.

There are no user-serviceable parts within. The compressor requires no lubrication or maintenance.

Device symbols

ON

OFF

Read the instructions

Type B

According to
93 / 42 / EEC Directive

Alternating current

Class II

The *Pro* nebulizer is suitable for use with all solutions and suspensions permitted for inhalation.

NEBULIZED SOLUTION

Budesonid

Albuterol

Beclametason Dipropionat

Tiopronin

Acetilcistein

Tiamfenicol Glicinat Acetilcistein

Flunisolid

Salbutamol

Ambroxol Clorid

Technical Data

Model:	MicroDrop® Pro
Electricity supply:	230 Volt, 50 HZ, 1A
Nebulization rate:	0.55 ml / min.
Particle size:	MMAD 1.8 µm (NaCl 0.9%) – API AEROSIZER MMAD 4.17 µm INSITEC – MALVERN
Nebulizer air flow:	3.3 l / min.
Compressor air flow:	7.5 l / min.
Compressor max. pressure:	2.50 bar
Nebulizer capacity:	max. 15 ml
Weight:	1.91 kg
Size:	(H) 106 x (L) 198 x (P) 233 mm
Noise level:	max. 52 dBA
Fuse:	T2.5 A, 250 V internal

Class II Device (Double Safety Insulation)

Type B Device (device with specific protection against electrical hazards)

Device for continuous use.

Device not protected against sprinkles.

Device not suitable for use in presence of anaesthetic substances inflammable with air, oxygen or nitrogen protoxide.

Device not suitable for use in anaesthesia or pulmonary ventilation systems.

Reference to standards: Electric Safety Standard CEI EN 60601-1
Electromagnetic Compatibility according to CEI EN 60601-1-2
Class IIa Medical Device according to 93 / 42 / EEC MDD

Manufacturer: Medel S.p.A., Via Micheli 7,
43056 S. Polo di Torrile (PR) Italy

Technical modifications reserved!

Ambient Conditions

Storing Conditions

Temperature: min. -25°C to max. 70°C
Humidity: min. 10% RH to max. 95% RH

Operating Conditions

Temperature: min. 10°C to max. 40°C
Humidity: min. 10% RH to max. 95% RH
Altitude: from 0 to 2000 m above sea level

Warranty certificate

Warranty conditions

The device is covered by a 5-year warranty from date of purchase for any defects in the design or in the materials used.

The warranty consists in the replacement and / or repair free of charge of originally defective components.

The warranty does not cover the accessories supplied with the device and those parts subject to normal wear and tear.

The device must be repaired by authorized technical service centres only.

The device must be sent to the authorized service centre for repair within 8 days after the defect is noticed.

Transport costs for the device are at user's charge.

Any repair out of warranty conditions are at user's charge.

The warranty does not cover failures resulting from misuse, abuse of the product or in case the damage is not due to manufacturer's liability (accidental fall, rough transportation, etc.).

The warranty does not involve any compensation for any kind of damages to persons or things, either direct or indirect, occurred when the device is not operating.

The warranty is valid from the date of purchase certified by the ticket or the invoice which have to be peremptorily attached to the warranty coupon.

The lack of warranty coupon properly filled in and validated by a certification of purchase voids the warranty.

Coupon to be returned in case of repairs

The warranty is valid only if accompanied by the ticket / invoice

Product-type:

Model:

Serial No.:

Date of purchase:

Buyer Data:

Name:

Address:

Phone number:

Defect description:

.....

.....

Signature for acceptance of the above mentioned warranty conditions

.....

Order Details

Description	REF
MicroDrop® Pro Inhalation device for youths and adults complete with compressor, <i>Pro</i> nebulizer for inhalation amounts up to 15 ml, mouth piece, nose piece, pressure tubing, mask adult, mask child, carrying bag and instruction manual	M 51803-00
Pro nebulizer Hospital Kit <i>Pro</i> nebulizer with mouth piece, mask adult, pressure tubing 2 m with fits-all connection, wall mount, clamp for wall rail	M 53510-01
Accessories	
Pro nebulizer incl. mouth piece, pressure tubing and instruction manual	M 53610-00
Nebulizer Set Pro comprising: <i>Pro</i> nebulizer, mouth piece with valve, nose piece, pressure tubing, mask adult, instruction manual	M 53610-01
Year Set Pro nebulizer comprising: <i>Pro</i> nebulizer with mouth piece, pressure tubing, Porex filter (air filter) for MicroDrop® compressor	M 53710-04
Mouth piece Pro nebulizer	M 50020-35200
Nose piece Pro nebulizer	M 50020-18500
Pressure tubing Pro nebulizer	M 50020-17900
Mask child Pro nebulizer	M 50020-32100
Mask adult Pro nebulizer	M 50020-32000
MicroDrop® Pro with mains cable	M 50010-38400
Porex filter air filter for MicroDrop® compressor	M 50010-30200
Carrying bag MicroDrop® Pro	M 50020-34700
Instruction manual MicroDrop® Pro compressor , multilingual	M 50020-35400
Instruction manual Pro nebulizer , multilingual	M 50020-35300

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This device fulfils the provisions of the EC Directive 93 / 42 / EEC.

Hersteller / Manufacturer (93 / 42 / EEC):
Designed, entwickelt und importiert von:
Designed, engineered and imported by:
MEDEL S.p.A.
S. Polo di Torrile (PR) Italy

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Kuhn und Bieri AG