



SALVACEL metal ring limb holder or immobilization system for non-agitated patients. The padded and breathable limb straps provide greater comfort and reduce pressure. Quick and easy installation with Velcro® closures. The metal rings provide additional adjustment for the wrists and ankles.

Made of soft-touch padded material (274gr/m²). 100% viscose inner coating, breathable 8mm polyester foam, and outer layer of smooth light blue velour fabric 100% polyamide, adhering to Velcro®. Polyester restraint straps with a breaking strength of up to 160Kg - 1560N. Quick closure system with Velcro®. Metal D-ring welded for durability.



LIMB HOLDER SYSTEM

SPECIFICATIONS

PRODUCT FAMILY: SALVACEL

PRODUCT NAME: Limb holder system

INTENDED USE: Preventive soft limb holder for non-agitated patients

BASIC UDI-DI: 843603535SVCN6

SIZE CHART

Ref.	Size	Measures (cm)	Straps length
SVC2476	M	30 x 8	108
SVC2476XL	L	35 x 8	108
SVC2476XL2	XL	35 x 8	119 (x2)

COLOUR

- Light blue
- Natural white

COMPOSITION

65% PES - 30% PA - 5% Viscose
Polyester straps
Nickel-free D shaped rings
LATEX FREE


LAUNDRY INSTRUCTIONS



- Machine wash 40°C
- Do not bleach
- Do not tumble dry
- Do not iron
- Do not dry clean

PACKAGING

- 1 Pair - Salvacel - Limb holder system



ATTENTION

- Single-patient product