

Foley Catheter Typ Tiemann

REF: 23.000.20.210 – 20.224

Unit: 10 pcs.

Packing Unit: 100 pcs.

Description

Catheters are used for controlled bladder release and postoperative drainage.

REF	Size (Fr/CH)	length	colour code	volume
23.000.20.210	10	40 cm	black	3 – 5 ml
23.000.20.212	12	40 cm	white	5 – 10 ml
23.000.20.214	14	40 cm	green	5 – 10 ml
23.000.20.216	16	40 cm	orange	5 – 10 ml
23.000.20.218	18	40 cm	red	5 – 10 ml
23.000.20.220	20	40 cm	yellow	5 – 10 ml
23.000.20.222	22	40 cm	purple	5 – 10 ml
23.000.20.224	24	40 cm	blue	5 – 10 ml

Material Composition

Product: Silicone, ABS, barium sulfate

Packing: individual packing: PE Bag, PET/PE, paper
box: cardboard
outer case: corrugated cardboard

Sterilisation

EO-Gas

23.000.20.212-23.000.20.224
QM/Rev.-Nr. 00/06.07.2010
QM/Rev.-Nr. 01/14.03.2013
QM/Rev. Nr. 02/28.01.2014
QM/Rev.-Nr. 03/28.04.2017
QM/Rev.-Nr. 04/23.09.2019
QM/Rev.-Nr. 05/26.09.2019