



SALVAFIX Memory PRO Complete system. The complete set is the global solution of secure restraint for bedridden patients. It permits having the maximum level of security of the patient and the medical staff, because of bed restraint is a not easy issue. As consists of individual elements can be used as required at the time.

All grouped in a sturdy carrying bag, capable of holding the entire system. It features a thermoplastic handle for a better grip and a closure with white injected zipper. It also includes tags for each size to be placed on the side of the bag, indicating the size of its contents.

COMPLETE SYSTEM

COMPONENTS



Harness



Wristbands



Sewn perineal band



Transportation bag



Waist Belt



Ankle bands



**Safety locks
Magnetic keys**



Colour code for sizes
S (Yellow) - M (Blue)
L (Green) - XL (Red)

SPECIFICATIONS

PRODUCT FAMILY: SALVAFIX Memory Pro

PRODUCT NAME: Complete system

INTENDED USE: Intended to be used as preventing patient self-injury and falls from the bed or stretcher

Basic UDI-DI: 843603535SVFNC

SIZE CHART

Ref.	Size
SVF3151P	S
SVF3150P	M
SVF3152P	L
SVF3153P	XL

COLOUR

- Natural white
- Light blue
- Brown

COMPOSITION

Polyester
Stainless steel eyelets
Stainless steel rings
LATEX FREE

LAUNDRY INSTRUCTIONS



- Machine wash 90°C
- Do not bleach
- Tumble dry low
- Do not iron

PACKAGING

Grouped in a carrying bag and fitted into a cardboard box:

- 1 harness with thorax band
- 1 bed strap for harness
- 1 waist belt
- 2 wristbands
- 2 ankle bands
- 1 bed strap for ankle bands
- 15 locks (16 for L and XL)
- 4 magnetic keys
- 1 carry bag
- 1 extendable key-ring
- Application stickers
- Instructions manual